

OCR Index As A Rapid Potency
Metric of Isolated Islets of Langerhans
By Chris Fraker
University of Miami School of Medicine

OBS: How it works


An oxygen sensitive Ruthenium-based fluorophore is adsorbed to silica particles and embedded in a silicone rubber matrix

- Diffusion from air into well and media

- Consumption by cells


- Diffusion across matrix

- O_2 - quenched $Ru(DPP)_3Cl_2$



METHODS:

Oxygen concentration is calculated using modified equations of Stern-Volmer Theory



t6

t3

$$I_0 / I_A = 1 + K_{SV} [O_2]_A$$

METHODS:

Oxygen concentration is calculated using modified equations of Stern-Volmer Theory



Helping all people
live healthy lives

- I_0 / I_A is the dynamic range of the plate fluorescence signal, DR
- K_{SV} is the Stern Volmer constant defined as $(DR-1) / [O_2]_A$
- Fluorescence values for each well were normalized (by division) to their initial value and then to the average value of no-cell controls at each time point (NRF)
- Ratio of DR/NRF = I_0 / I_{NRF}
- $DR/NRF = 1 + K_{SV} * [O_2]_{NRF}$
- So, $NRF = DR / (1 + K_{SV} * [O_2]_{NRF})$
- Fluorescence values so-normalized were used to compute oxygen concentrations per Equation [1].

$$[O_2]_{NRF} = \frac{DR}{NRF-1} \times \frac{1}{K_{SV}} \quad [1]$$

t3

Transient State dO_2 Slope and Equilibrium OCR Calculations



Short-term (< ~1 hr)

- Quantitate based on rate of change of signal
- Transient dO_2 Slope ~ Slope of $1-1/NRF$ vs. time
- Ideal for Islets

Steady State

- Quantitated based on equilibrium OBS signal
- $OCR = D S L \Delta p / h$
- $OCR \sim Actual \ 1-1/NRF \sim DdO_2$
- Additional comparative metric of islet potency. "True OCR"
- Time to equilibrium depends on volume of medium used.

$$T = h^2/D$$

$$D = 3.3 \times 10^{-5} \text{ cm}^2/\text{s}$$

$$S = 0.31 \text{ cm}^2 \text{ for a 96-well plate}$$

$$h = 0.65 \text{ cm for } 200 \mu\text{L}$$

$$L = 6.0 \times 10^7 \text{ L/cm}^3 * \text{fmol}/\mu\text{mol} * \text{s/min}$$




$$\Delta p = \text{in } \mu\text{M.}$$

METHODS:

Triplicate wells of IEQs in Low (2.2mM) and High (22mM) Glucose Supplemented Krebs Buffer



Helping all people live healthy lives



Islets Low glucose

Islets High glucose

Media CNT

Sulfite CNT


- Measure fluorescence change for 200 cycles 5 minute intervals
- Calculate transient steady state t_{60} dO_2 and equilibrium OCR and normalize for DNA content.
- Index = Mean High glucose values/ Mean Low glucose values



t10 Example Assessment: Islet Shipment From UPS/Seattle 6-29-05




Helping all people
live healthy lives



Transient Period Index: 1.51

Equilibrium Index: 1.56



Transient period dO2 slope and Equilibrium OCR scale kinetically with cell number



dO₂ and Equilibrium OCR Scale Linearly with IEQ



Helping all people live healthy lives




Both metrics of OCR scale with IEQ and each other




Helping all people live healthy lives

t9

Rapid dO₂/DNA vs. Equilibrium OCR/DNA



Rapid dO₂ raw vs. Equilibrium OCR raw



Diabetic Mouse Model for Potency



Helping all people
live healthy lives


- Chemically Induced Diabetic Mice
- Islets Transplanted into Kidney Capsule
- Recipients chosen based on body weight and blood glucose profile similarities.
- 2000 IEQ per mouse
- Monitor blood glucose^{t11} daily
- Cure is defined as two consecutive days < 200 mg/dl fasting blood glucose
- Days to cure used as comparative measure of^{t12} potency.

OCR/DNA shows some correlation to reversal...




Helping all people
live healthy lives

t13



...but OCR Index is much more predictive




- ◆ solid red = pelleted --> all failed
- ◇ open red = not pelleted, but failed at least some replicates
- ◆ yellow = 100% reversal, but not all by day 3
- ◆ green = 100% reversal by day 3
- ◇ open green: n=1
- * asterisk = PFC

Kaplan-Meier type survival analysis.

- P-values calculated using a log-rank test. $P < .00001$
- Top curve "cure" curve for mice with ratios > 1.27
- Bottom curve, mice with ratios < 1.27 .
- Dashed lines 95% confidence intervals for the fraction cured.



Helping all people
live healthy lives



In 31 human islet preparations analyzed, statistical analysis bears out that index is by far the best predictor

t14



Acknowledgements:

BD

- Mark Timmins, PhD.
- Richard D. Guarino, MS.
- Sharon Presnell, PhD.
- Perry D Haaland, PhD.

University of Minnesota Diabetes Institute for Immunology and Transplant

- Klearchos Papas, PhD.

University of Miami Diabetes Research Institute

- Damaris Molano, DVM
- Antonello Pileggi, MD
- Hirohito Ichii, MD, PhD.
- Luca Inverardi, MD
- Rafaella Poggioli, MD
- Elsie Zahr
- Jill Collins
- Camillo Ricordi, MD



Helping all people
live healthy lives